ULTRA-PRO 3™

NEEDLE GUIDANCE SYSTEM


REFERENCE GUIDE

Lock
Zaključavanje
Zámek
Lås
Vergrendeling
Lukk
Lukitus
Verrou
Riegel
Κλειδωμα
נעילה
Zár
Kunci
Fermo
ロック
잠금

Insert
Umetak
Vložka
Indsats
Insteekhuls
Sisestusosa
Sisäkäppäle
Insert
Einsatz
Εισαγωγή
תוספת
Bevezetõ
Masukkan
Inserto
インサート
인서트


Turquoise Indicator
Tirkizni indikator
Tyrkysová značka
Turkis indikator
Turquoise indicator
Türkiissinine indikaator
Turkoosi merkki
Indicateur turquoise
Türkisfarbener Indikator
Ένδειξη χρώματος
πικουάζ
מחון טורקיז
Türkizkék indikátor
Indikator Pirus
Indicatore turchese
青緑色インジケータ
청록색 표시기

Needle Guide Body
Tijelo Vodilice Igle
Tělo vodiče jehly
Nåleguidens krop
Naaldvoerder
Nõela juhikere
Neulanohjaimen runko
Corps d'aiguille-guide
Nadelführungskörper
Σώμα οδηγού βελόνας
גוף מוביל מחט
Tűvezetõ
Isi Pengarah Jarum
Corpo della guida dell'ago
ニードルガイド本体
니들 가이드 본체

Turquoise Quick-Release Lever
Tirkizna poluga za brzo otpuštanje
Tyrkysová rychloupínací páčka
Turkis hurtigudløsningsstap
Turquoise snelontkoppelingshefboom
Türkiissinine kiirvabastuse hoob
Turkoosi pikavapautusvipu
Lever à blocage rapide turquoise
Türkisfarbener Schnelltrennhebel
Μοχλός ταχείας απελευθέρωσης χρώματος
πικουάζ
מנוף שחרור מהיר בצבע טורקיז
Türkizkék gyorskiengedő kar
Tuas Pelepas Cepat Pirus
Leva turchese del dispositivo di sblocco rapido
青緑色クイックリリースレバー
청록색 신속 배출 레버

Fiksators	Ieliktnis	Tirkīzkrāsas indikators
Fiksatorius	[dedamoji dalis	Turkio rodiklis
Lāš	Innlegg	Turkisfarget indikator
Blokada	Wkładka	Wskaźnik turkusowy
Trava	Inserir	Indicador turquesa
Fecho	Insertores	Indicador de cor turquesa
Blocaj	Insertie	Indicator turcoaz
Фиксатор	Вставка	Бирюзовый индикатор
锁销	插件	蓝绿色指示器
Poistka	Vložka	Tyrkysový indikátor
Cierre	Inserción	Indicador turquesa
Lāš	Införare	Turkosfärgad indikator
ตัวล็อก	ส่วปลาย	ไฟแสดงสถานะสี เทอร์
鎖	插件	คอยซ์
Kilit	Ek Parça	藍綠色指示器
Khóa	Chi tiết đệm	Turkuaz Gösterge
		Chỉ Báo Màu Ngọc Lam

Adatas vadīklas korpus
 Adatos kreipiamoji dalis
 Nåleføringsenhet
 Korpus prowadnicy igły
 Corpo da guia da agulha
 Corpo da guia para agulhas
 Corp ghidaj ac
 Корпус проводника иглы
 针导引器体
 Telo vodidla ihly
 Guia para aguja
 Nålguidehölje
 ตัววอ ปรกรณ์ ควค มเช้ มเจาะ
 針導引器主體
 İğne Kılavuzu Gövdesi
 Thân Thanh Dẫn Kim Tiêm

Tirkīzkrāsas ātrās atbrīvošanas svira
 Turkio greito atlaisvinimo svirtelė
 Turkisfarget hurtigutløsningsspak
 Turkusowa dźwignia szybkiego zwalniania
 Alavanca de liberação rápida turquesa
 Alavanca de libertação rápida de cor turquesa
 Mâner de eliberare rapidă turcoaz
 Бирюзовый рычаг быстрого высвобождения
 蓝绿色快速释放杆
 Tyrkysový rýchlopínací držiak
 Pestaña turquesa de liberación rápida
 Turkosfärgad snabbutlösningsspak
 ตัว โยกปลดเร็ว วสี เทอร์ คอยซ์
 藍綠色快速釋放桿
 Turkuaz Hızlı Serbest Birakma Kolu
 Cần Nhà Nhanh Màu Ngọc Lam

Symbol	Title of Symbol	Description of Symbol
	Manufacturer (ISO 15223-1, 5.1.1)	Indicates the medical device manufacturer.
	Authorized Representative in the European Community (ISO 15223-1, 5.1.2)	Indicates the Authorized Representative in the European Community.
	Date of manufacture (ISO 15223-1, 5.1.3)	Indicates the date when the medical device was manufactured.
	Use-by date (ISO 15223-1, 5.1.4)	Indicates the date after which the medical device is not to be used.
	Batch code (ISO 15223-1, 5.1.5)	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog number (ISO 15223-1, 5.1.6)	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Sterilized using ethylene oxide (ISO 15223-1, 5.2.3)	Indicates a medical device that has been sterilized using ethylene oxide.
	Do not use if package is damaged (ISO 15223-1, 5.2.8)	Indicates a medical device that should not be used if the package has been damaged or opened.
	Keep away from sunlight (ISO 15223-1, 5.3.2)	Indicates a medical device that needs protection from light sources.
	Do not reuse (ISO 15223-1, 5.4.2)	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Consult instructions for use (ISO 15223-1, 5.4.3)	Indicates the need for the user to consult the instructions for use.
	Not made with natural rubber latex (ISO 15223-1, 5.4.5 and Annex B)	Indicates that natural rubber latex was not used in the manufacturing of the product, its container, or its packaging.
	European Conformity (EU MDR 2017/745, Article 20)	Indicates manufacturer declaration that the product complies with the essential requirements of the relevant European health, safety and environment protection legislation.
	Medical Device (MedTech Europe Guidance: Use of Symbols to Indicate Compliance with the MDR)	Indicates the product is a medical device.
	Quantity (IEC 60878, 2794)	To indicate the number of pieces in the package.

BRACKET / NEEDLE GUIDE**INTENDED USE**

The bracket and disposable guide provide physicians with a tool for performing needle-guided (or catheter) procedures with the use of diagnostic ultrasound transducers.

COVER**INTENDED USE**

Protective cover or sheath placed over diagnostic ultrasound transducer / probe / scanhead instruments. The cover allows use of the transducer in scanning and needle guided procedures for body surface, endocavity, and intra-operative diagnostic ultrasound, while helping to prevent transfer of microorganisms, body fluids, and particulate material to the patient and healthcare worker during reuse of the transducer (both sterile and non-sterile covers). The cover also provides a means for maintenance of a sterile field (sterile covers only). CIVCO Poly Ultrasound Transducer Covers are furnished sterile & non-sterile; single use patient / procedure, disposable.

INDICATIONS FOR USE

- Abdominal - Biopsy, Fine-Needle Aspiration, Regional Anesthesia Injections.
- Regional Anesthesia - Nerve Block, Musculoskeletal Injections.
- Small Parts - Biopsy, Fine-Needle Aspiration, Regional Anesthesia Injections.

PATIENT POPULATION

General Purpose Guidance Systems are intended for use with appropriately sized transducer intended for use on adults of all body habitus, including small or slight frame.

INTENDED USERS

General Purpose Guidance Systems should be used by clinicians medically trained in the use and interpretation of ultrasound for interventions. There are no other unique skills or user abilities required for Bracket and Needle Guide use. This may include, but is not limited to: Anesthesiologists, Interventional Radiologists, Nephrologists, Radiologists, Sonographers and Surgeons.

PERFORMANCE CHARACTERISTICS

- Ultra-Pro 3™ Needle Guidance System provides accurate needle path with on-screen guidelines, supports wide range of gauge sizes, and incorporates quick-release feature for detachment of guide.
- Cover serves as a viral barrier to protect patients, users, and equipment from cross-contamination.

NOTE: For a summary of clinical benefits for this product, visit www.CIVCO.com.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

- Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's user guide.
- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.
- See your system's user guide for reprocessing transducer between use.
- Review product labeling for sterile symbol. If sterile symbol appears on label, cover and guide are sterilized using ethylene oxide.
- Sterile product is intended for single-use only.
- Do not use if integrity of packaging is violated.
- Do not use if expiration date has passed.
- Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.
- CIVCO's needle guidance systems are designed and validated as an integrated two-part system consisting of a custom reusable bracket and a disposable guide. Use of non-CIVCO guides on a CIVCO bracket or locating feature may result in serious injury to the patient and/or user. To ensure safe and effective use in compliance with system software guidelines, only CIVCO needle guides are recommended.
- Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.

- For illustration purposes only, transducer and bracket may be shown without a transducer cover. Always place a cover over transducer and bracket to protect patients and users from cross-contamination.
- Do not use if needle guide is damaged or does not fit properly.
- Ensure no air bubbles are present prior to use. Air left between cover and transducer lens may cause poor image quality.
- If the product malfunctions during use or is no longer able to achieve its intended use, stop using the product and call CIVCO.
- Report serious incidents related to the product to CIVCO and the competent authority in your Member State or appropriate regulatory authorities.

NOTE: Product is not made with natural rubber latex.


NEEDLE PATH VERIFICATION

- Before using component for first time, perform needle path verification to verify system and biopsy guide relationships, as described in system user guide.
- Insert may be adjusted once needle is in guide to achieve desired feel.

NOTE: Use appropriate needle length to reach target area.


USING NEEDLE GUIDANCE SYSTEM

1. Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
2. Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
3. Secure cover with enclosed bands.
4. Inspect cover to ensure there are no holes or tears.


SELECT INSERT

1. Select appropriate size needle guide insert based on gauge size required.


ATTACHING INSERT ONTO NEEDLE GUIDE


1. Slide needle gauge insert onto turquoise lever.

2. Align gauge size with turquoise indicator.


3. Snap turquoise quick-release lever closed.


ATTACH AND LOCK NEEDLE GUIDE TO BRACKET

1. Using proper sterile technique, snap unlocked needle guide onto attachment area of bracket.
2. Push lock into locked position.
3. Verify gauge selection.

ACTIVATING NEEDLE GUIDE QUICK-RELEASE

1. Press turquoise lever towards bracket to activate needle guide quick-release function.

STORAGE CONDITIONS

- Avoid storing product in areas of temperature extremes or in direct sunlight.
- Store in a cool, dry place.

DISPOSAL

WARNING

- *Dispose of single-use components as infectious waste. Clean and disinfect reusable components after each use.*

NOTE: For questions or to order additional CIVCO products, please call +1 319-248-6757 or 1-800-445-6741 or visit www.CIVCO.com.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ULTRA-PRO 3 IS A TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.
MPS GmbH
RQ Solutions Pty Ltd
CIVCO Medical Instruments Co., Inc.