OMNITRAX™

ACTIVE PATIENT TRACKER


For use in MR environment

REFERENCE GUIDE

For use with
Canon Aplio i700, i800, i900;
Esaote MyLab 9, MyLab Class, MyLab Eight, MyLab Twice;
FUJIFILM ARIETTA 850;
GE LOGIQ E9 XDclear 2.0, LOGIQ E10;
SuperSonic Imagine Aixplorer® version 11.2 or later
transducers

CE
1639

Symbol	Title of Symbol	Description of Symbol
	Manufacturer (ISO 15223-1, 5.1.1)	Indicates the medical device manufacturer.
	Authorized Representative in the European Community (ISO 15223-1, 5.1.2)	Indicates the Authorized Representative in the European Community.
	Date of manufacture (ISO 15223-1, 5.1.3)	Indicates the date when the medical device was manufactured.
	Use-by date (ISO 15223-1, 5.1.4)	Indicates the date after which the medical device is not to be used.
	Batch code (ISO 15223-1, 5.1.5)	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog number (ISO 15223-1, 5.1.6)	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	MR Safe (ASTM F2503, 6.4.4)	Indicates device is safe for use in MR environment.
	Do not use if package is damaged (ISO 15223-1, 5.2.8)	Indicates a medical device that should not be used if the package has been damaged or opened.
	Keep away from sunlight (ISO 15223-1, 5.3.2)	Indicates a medical device that needs protection from light sources.
	Single patient multiple use (ISO 7000, 3706)	Indicates the medical device may be used multiple times (multiple procedures) on a single patient.
	Consult instructions for use (ISO 15223-1, 5.4.3)	Indicates the need for the user to consult the instructions for use.
	Not made with natural rubber latex (ISO 15223-1, 5.4.5 and Annex B)	Indicates that natural rubber latex was not used in the manufacturing of the product, its container, or its packaging.
	European Conformity (EU MDR 2017/745, Article 20)	Indicates manufacturer declaration that the product complies with the essential requirements of the relevant European health, safety and environment protection legislation.
	Medical Device (MedTech Europe Guidance: Use of Symbols to Indicate Compliance with the MDR)	Indicates the product is a medical device.
	Quantity (IEC 60878, 2794)	To indicate the number of pieces in the package.

omniTRAX™ ACTIVE PATIENT TRACKER**INTENDED USE**

The device is intended to provide physicians with a tool for electromagnetic (EM) tracking of instruments with respect to pre-acquired or real-time data. The device is intended for use in Magnetic Resonance Imaging (MRI) and ultrasound (US) environments.

Device is MR Safe.

INDICATIONS FOR USE

- Abdominal - Image guided minimally invasive puncture procedures.
- Small Parts - Image guided minimally invasive puncture procedures.
- Transrectal - Image guided minimally invasive puncture procedures.

PATIENT POPULATION

omniTRAX™ Active Patient Tracker is intended for use on adults of all body habitus, including small or slight frame.

INTENDED USERS

omniTRAX™ Active Patient Tracker should be used by clinicians medically trained in the use and interpretation of ultrasound for interventions. There are no other unique skills or user abilities required for device use. This may include, but is not limited to: Interventional Radiologists, Radiologists, Sonographers and Surgeons.

PERFORMANCE CHARACTERISTICS

The omniTRAX™ Active Patient Tracker provides clinicians automatic image registration of fused images for uninterrupted navigation when using real-time ultrasound with previously acquired MR volume data sets.

NOTE: For a summary of clinical benefits for this product, visit www.CIVCO.com.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.


WARNING

- *Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's user guide.*
- *Prior to use, refer to your system's user guide for use of general purpose sensor.*
- *Prior to use, inspect device for signs of damage, if damage is evident do not use.*
- *Prior to use, ensure general purpose sensor is tracking properly.*
- *Active patient tracker is packaged non-sterile and is intended for single-patient use only.*
- *To avoid possible position distortion, keep general purpose sensor away from metal.*
- *Ensure device is securely attached, failure to do so could result in inaccuracy.*
- *Ensure device is secure on patient prior to use.*
- *Ensure general purpose sensor is securely attached to the active patient tracker.*
- *If device leakage occurs, fluid is non-cytotoxic and routine hospital clean-up protocol should be followed.*
- *If the product malfunctions during use or is no longer able to achieve its intended use, stop using the product and call CIVCO.*
- *Report serious incidents related to the product to CIVCO and the competent authority in your Member State or appropriate regulatory authorities.*


NOTE: Product is not made with natural rubber latex.

ATTACHING omniTRAX™ ACTIVE PATIENT TRACKER TO PATIENT


NOTE: MR image quality may be compromised if the anatomical region of interest is in close proximity of the device. Image artifact extends approximately 2mm from the adhesive pads of the device. To optimize imaging parameters, it is recommended to position the device away from anatomical region of interest, but within the range of the MR scan.

1. Remove adhesive liners and position active patient tracker onto patient's skin.
2. Mark device location on patient's skin using active patient tracker marking locations.

**ATTACHING GENERAL PURPOSE SENSOR TO ACTIVE PATIENT TRACKER**

1. Place general purpose sensor onto sensor mount area.
2. Rotate general purpose sensor into active patient tracker to secure.
3. Following proper sterile technique, place a sterile patient drape over patient, ensuring active patient tracker and general purpose sensor are positioned under drape.

**DISPOSAL****WARNING**

- When active patient tracker can no longer be securely attached to patient or becomes contaminated, dispose as infectious waste.

REPROCESSING**WARNING**

- Refer to general purpose sensor instruction for reprocessing.

STORAGE CONDITIONS

- Avoid storing product in areas of temperature extremes or in direct sunlight.
- Store in a cool, dry place.

NOTE: For questions or to order additional CIVCO products, please call +1 319-248-6757 or 1-800-445-6741 or visit www.CIVCO.com.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. OMNITRAX IS A TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.
MPS GmbH
RQ Solutions Pty Ltd
CIVCO Medical Instruments Co., Inc.