


Verza™

GUIDANCE SYSTEM

REFERENCE GUIDE

For use with
Esaote
LMX 4-20
transducers

Refer to separate insert for Symbol Glossary.

GUIDANCE SYSTEM

INTENDED USE

The system provides guidance for precise instrument placement of common interventional devices by positioning the device relative to the ultrasound transducer and the resulting image during a diagnostic or therapeutic procedure. This guidance system is intended for use with pediatric and adult patients.

COVER

INTENDED USE

Protective cover or sheath placed over diagnostic ultrasound transducer / probe / scanhead instruments. The cover allows use of the transducer in scanning and needle guided procedures for body surface, endocavity, and intra-operative diagnostic ultrasound, while helping to prevent transfer of microorganisms, body fluids, and particulate material to the patient and healthcare worker during reuse of the transducer (both sterile and non-sterile covers). The cover also provides a means for maintenance of a sterile field (sterile covers only). CIVCO Poly Ultrasound Transducer Covers are furnished sterile & non-sterile; single use patient / procedure, disposable.

INDICATIONS FOR USE

- Abdominal - Biopsy, Fine-Needle Aspiration, Regional Anesthesia Injections.
- Regional Anesthesia - Nerve Block, Musculoskeletal Injections.
- Small Parts - Biopsy, Fine-Needle Aspiration, Regional Anesthesia Injections.

PATIENT POPULATION

The Verza Guidance System is intended for use with pediatric and adult patients.

INTENDED USERS

General purpose guidance systems should be used by clinicians medically trained in the use and interpretation of ultrasound for interventions. There are no other unique skills or user abilities required for Bracket and Needle Guide use. This may include, but is not limited to: Anesthesiologists, Interventional Radiologists, Nephrologists, Radiologists, Sonographers and Surgeons.

PERFORMANCE CHARACTERISTICS

- The Verza Needle Guidance System's innovative design reduces blind zone and unguided needle length while providing high stability for improved accuracy. Verza offers maximum versatility with a five-angle approach for improved anatomical access while also featuring an expanded gauge range, accommodating 14 different sizes of interventional instruments. Verza features easy-to-read gauge sizes, angle indicator and repeatable bolt action quick-release for secure instrument detachment.
- The Verza Bracket design includes a simplified multi-angle design free of lock pins and removable parts. The bracket also offers a low-profile fit, sleek design for ease of cleaning, and secure user-friendly closure door for ease of application and removal.
- Cover serves as a viral barrier to protect patients, users, and equipment from cross contamination.

NOTE: For a summary of clinical benefits for this product, visit www.CIVCO.com.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

- Verza Guidance System is designed and validated for use with transducer(s) listed on the cover of this instruction.
- Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your ultrasound system's user guide.
- Review product labeling for sterile symbol. If sterile symbol appears on label, cover and guide are sterilized using ethylene oxide.
- Sterile product is intended for single-use only.
- Do not use if integrity of packaging is violated.
- Do not use if expiration date has passed.
- Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.
- CIVCO's guidance systems are designed and validated as an integrated two-part system consisting of a custom reusable bracket and a disposable guide. Use of non-CIVCO guides on a CIVCO bracket or locating feature may result in serious injury to the patient and/or user. To ensure safe and effective use in compliance with system software guidelines, only CIVCO guides are recommended.
- Bracket is packaged non-sterile and is reusable. To avoid possible patient contamination, ensure bracket is properly cleaned and disinfected before each use.
- Do not use if device appears damaged or does not fit properly.
- Do not use if locating feature is damaged or guide does not fit properly.
- Ensure gauge insert is in closed position prior to beginning procedure.
- Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.
- For illustration purposes only, transducer and bracket may be shown without a transducer cover. Always place a cover over transducer and bracket to protect patients and users from cross-contamination.
- Ensure no air bubbles are present prior to use. Air left between cover and transducer lens may cause poor image quality.
- If the product malfunctions during use or is no longer able to achieve its intended use, stop using the product and call CIVCO.
- Report serious incidents related to the product to CIVCO and the competent authority in your Member State or appropriate regulatory authorities.

NOTE: Product is not made with natural rubber latex.

NEEDLE PATH VERIFICATION


- Before using component for first time, perform needle path verification to verify system and Verza guide relationships, as described in ultrasound system user guide.
- Before each use, verify guidance system angles (5) correspond with system software.

NOTE: Use appropriate needle length to reach target area. Verza guide channel measures 3.4cm (1.34in).


USING GUIDANCE SYSTEM

1. With latch open, align bracket with locating features on transducer.

2. Close latch around transducer and snap to secure. Ensure bracket is securely attached.

3. Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
4. Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
5. Secure cover with enclosed bands.
6. Inspect cover to ensure there are no holes or tears.
7. Remove Verza from guide storage column of gauge insert tree. Using proper sterile technique, place Verza guide onto VerzaLink locating feature.


8. Rotate guide locking lever into locked position. Ensure guide is securely attached.
9. **FIXED ANGLE APPROACH:**
Using angle selector grip, slide angle indicator window to desired position and push angle selector lock toward transducer to secure.

NOTE: Ensure angle selector lock is secure in desired position.

FREE ANGLE APPROACH:


Leave angle selector lock in unlocked position and instrument remains in-plane within constraints of guide angle range.

10. Select appropriate size gauge insert and slide into guide body chute.


BOLT-ACTION QUICK-RELEASE

1. Rotate quick-release tab to open guide insert funnel.

REMOVING GUIDANCE SYSTEM

1. Rotate guide locking lever to unlocked position to remove guide from VerzaLink locating feature of bracket.
2. Lift bracket latch to remove from transducer. Retain bracket for reprocessing.

REPROCESSING

WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.
- Dispose of single-use components as infectious waste. Clean and disinfect reusable components after each use.
- See your ultrasound system's user guide for reprocessing transducer between use.
- This procedure for reprocessing has been validated for effectiveness and compatibility. Device could be damaged or cross-contaminated due to improper reprocessing.
- Do not gas sterilize or autoclave bracket.
- If bracket damage is evident, discontinue use and contact your CIVCO representative.

CLEANING

1. After each use, discard single-use, disposable components and keep bracket from drying out until complete cleaning can be accomplished.
2. Inspect bracket for damage such as cracks or breakage.

MANUAL CLEANING

NOTE: Manual cleaning efficacy of this process was proven using Enzo® Enzymatic Detergent.

1. Rinse bracket under cold running potable water for one minute to remove gross soil. During rinse, brush all surfaces with soft-bristled brush while actuating moveable parts. Ensure all cracks, crevices, and hard-to-reach areas are brushed.
2. Prepare detergent solution following manufacturer's instructions for use.
3. Fully immerse bracket in prepared detergent solution and allow to soak at least the minimum contact time listed in solution manufacturer's instructions for use. Brush all surfaces with soft-bristled brush while actuating moveable parts. Ensure all cracks, crevices, and hard-to-reach areas are brushed.
4. If visible soil still exists, repeat Step 3.
5. Rinse bracket under cold running potable water for one minute.
6. Dry with clean, lint-free cloth. Visually inspect to ensure all surfaces are clean and dry.

AUTOMATED CLEANING

NOTE: Automated cleaning efficacy of this process was proven using neodisher® Mediclean forte.

1. Rinse bracket under cold running potable water for one minute to remove gross soil. During rinse, brush all surfaces with soft-bristled brush while actuating moveable parts. Ensure all cracks, crevices, and hard-to-reach areas are brushed.
2. Prepare detergent solution following manufacturer's instructions for use.
3. Fully immerse bracket in prepared detergent solution and allow to soak at least the minimum contact time listed in solution manufacturer's instructions for use. During rinse, brush all surfaces with soft-bristled brush while actuating moveable parts. Ensure all cracks, crevices, and hard-to-reach areas are brushed.
4. If visible soil still exists, repeat Step 3.
5. Place bracket onto rack system in washer, following hospital protocol for washer load requirements.
6. Proceed with appropriate phase cycles using parameters listed below:

Phase	Recirculation Time (Minutes)	Temperature	Detergent Type and Concentration (if applicable)
Pre-wash 1	02:00	Cold potable water	N/A
Wash 1	01:00	43°C potable water (set point)	neodisher® Mediclean forte - 9.9ML/L
Rinse 1	01:00	43°C potable water (set point)	N/A
Dry Time	07:00	90°C	N/A

7. Remove bracket from washer. Visually inspect to ensure all surfaces are clean and dry.

DISINFECTION

LOW-LEVEL DISINFECTION - SANI-CLOTH® AF3 GERMICIDAL DISPOSABLE WIPE

NOTE: Low-Level disinfection efficacy of this process was proven using Sani-Cloth® AF3 Germicidal Disposable Wipe.

1. Wipe bracket with Sani-Cloth® AF3 Wipe thoroughly to remove gross soil.

2. With new Sani-Cloth® AF3 Wipe, thoroughly wet surface for at least the minimum time listed in disinfectant manufacturer's instructions for use. Use additional wipes as needed. Focus on cracks, crevices, and hard-to-reach areas, including latch.
3. Thoroughly rinse bracket under running potable water for two minutes. While rinsing, wipe with clean, lint-free cloth to remove disinfectant residuals.
4. Allow bracket to dry. Visually inspect to ensure all surfaces are clean and dry.

HIGH-LEVEL DISINFECTION - REVITAL-OX® RESERT® HIGH LEVEL DISINFECTANT

NOTE: High-Level disinfection efficacy of this process was proven using Revital-Ox® RESERT® High Level Disinfectant.

1. Prepare high-level disinfectant according to manufacturer's instructions.
2. Fully immerse bracket in disinfectant solution and agitate to remove all air bubbles from surface.
3. Allow bracket to soak in disinfectant solution for at least the minimum contact time listed in disinfectant manufacturer's instructions for use.
4. Rinse thoroughly by fully immersing bracket in large volume of potable water, agitate, and allow to set for at least one minute.
5. Repeat Step 4 with new bath of potable water.
6. Remove bracket from potable water.
7. Dry with sterile, lint-free cloth. Visually inspect to ensure all surfaces are clean and dry.

HIGH-LEVEL DISINFECTION - CIDEX® OPA SOLUTION

NOTE: High-Level disinfection efficacy of this process was proven using CIDEX® OPA Solution.

1. Prepare high-level disinfectant according to manufacturer's instructions.
2. Fully immerse bracket in disinfectant solution. Use syringe to flush and remove air bubbles. Focus on cracks, crevices, and hard-to-reach areas. Actuate all moveable parts.
3. Allow bracket to soak in disinfectant solution for at least the minimum contact time listed in disinfectant manufacturer's instructions for use.
4. Rinse thoroughly by fully immersing bracket in large volume of potable water, agitate, and allow to set for at least one minute.
5. Remove bracket from potable water.
6. Dry with sterile, lint-free cloth. Visually inspect to ensure all surfaces are clean and dry.

HIGH-LEVEL DISINFECTION - GIGASEPT® FF (NEW)

NOTE: High-Level disinfection efficacy of this process was proven using gigasept® FF (new).




1. Prepare high-level disinfectant according to manufacturer's instructions.
2. Fully immerse bracket in disinfectant solution. Use syringe to flush and remove air bubbles. Focus on cracks, crevices, and hard-to-reach areas. Actuate all moveable parts.
3. Allow bracket to soak in disinfectant solution for at least the minimum contact time listed in disinfectant manufacturer's instructions for use.
4. Thoroughly rinse bracket under running potable water for one minute. Remove air bubbles by agitating and flushing with syringe.
5. Dry with sterile, lint-free cloth. Visually inspect to ensure all surfaces are clean and dry.

STORAGE CONDITIONS

- Avoid storing product in areas of temperature extremes or in direct sunlight.
- Store in a cool, dry place.

NOTE: Bracket is validated for 5-year product life.

NOTE: For questions or to order additional CIVCO products, please call +1 319-248-6757 or 1-800-445-6741 or visit www.CIVCO.com.

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COPYRIGHT © 2023 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. VERZA AND VERZALINK ARE TRADEMARKS OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MPS GmbH
CIVCO Medical Instruments Co., Inc.