ENVISION™

100% Gel-Free Ultrasound ProceduresSM

REFERENCE GUIDE

CE
1639

Symbol	Title of Symbol	Description of Symbol
	Manufacturer (ISO 15223-1, 5.1.1)	Indicates the medical device manufacturer.
	Authorized Representative in the European Community (ISO 15223-1, 5.1.2)	Indicates the Authorized Representative in the European Community.
	Date of manufacture (ISO 15223-1, 5.1.3)	Indicates the date when the medical device was manufactured.
	Use-by date (ISO 15223-1, 5.1.4)	Indicates the date after which the medical device is not to be used.
	Batch code (ISO 15223-1, 5.1.5)	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog number (ISO 15223-1, 5.1.6)	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Sterilized using ethylene oxide (ISO 15223-1, 5.2.3)	Indicates a medical device that has been sterilized using ethylene oxide.
	Do not use if package is damaged (ISO 15223-1, 5.2.8)	Indicates a medical device that should not be used if the package has been damaged or opened.
	Keep away from sunlight (ISO 15223-1, 5.3.2)	Indicates a medical device that needs protection from light sources.
	Do not reuse (ISO 15223-1, 5.4.2)	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Consult instructions for use (ISO 15223-1, 5.4.3)	Indicates the need for the user to consult the instructions for use.
	Not made with natural rubber latex (ISO 15223-1, 5.4.5 and Annex B)	Indicates that natural rubber latex was not used in the manufacturing of the product, its container, or its packaging.
	European Conformity (EU MDR 2017/745, Article 20)	Indicates manufacturer declaration that the product complies with the essential requirements of the relevant European health, safety and environment protection legislation.
	Medical Device (MedTech Europe Guidance: Use of Symbols to Indicate Compliance with the MDR)	Indicates the product is a medical device.
	Quantity (IEC 60878, 2794)	To indicate the number of pieces in the package.

INTENDED USE

The Envision cover and Envision pad are intended as a microbial barrier between the patient and medical imaging electronics.

INDICATIONS FOR USE

- Abdominal - Diagnostic imaging and minimally invasive puncture procedures.
- Small Parts - Diagnostic imaging and minimally invasive puncture procedures.
- Regional Anesthesia - Minimally invasive puncture procedures.
- Vascular Access - Vessel identification and catheter placement.
- Surgical - Diagnostic imaging and puncture procedures.

PATIENT POPULATION

Envision cover and Envision pad are intended for use with all body size and habitus.

INTENDED USERS

Transducer covers should be used by clinicians medically trained in the proper methods of sterile technique. There are no other unique skills or user abilities required for Transducer Cover use. User groups may include, but are not limited to: Anesthesiologists, Department Nursing staff, Emergency Medicine physicians, Endocrinologists, Interventional Radiologists and Interventional Radiology Technologists, Radiologists, Sonographers, Surgeons, Surgical Technicians, and other clinicians trained in the use of ultrasound.

PERFORMANCE CHARACTERISTICS

Envision cover and Envision pad serve as a viral barrier and enable 100% gel-free ultrasound procedures.

NOTE: For a summary of clinical benefits for this product, visit www.CIVCO.com.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

- Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's user guide.
- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.
- See your system's user guide for reprocessing transducer between use.
- Product is intended for single-use only.
- Do not use if integrity of packaging is violated.
- Do not use if expiration date has passed.
- If the product malfunctions during use or is no longer able to achieve its intended use, stop using the product and call CIVCO.
- Report serious incidents related to the product to CIVCO and the competent authority in your Member State or appropriate regulatory authorities.
- Review product labeling for sterile symbol. If sterile symbol appears on label, cover is sterilized using ethylene oxide.
- Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.
- Always place a cover over transducer to protect patients and user from cross-contamination.
- Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.
- Organic solvents, such as alcohol, may damage the coating.
- Not intended for use with the central circulatory system or central nervous system including the dura.

NOTE: Product is not made with natural rubber latex.

NOTE: THIS PRODUCT DOES **NOT** REQUIRE THE USE OF ULTRASOUND GEL.

APPLYING THE COVER

NOTE: Instruction is for a two-person cover application.

1. **Sterile user:** Insert hand completely into cover. CIVCO liner and images should be facing upward.



2. **Sterile user:** Peel liner away from cover to expose adhesive.

NOTE: If part of cover sticks to adhesive, carefully peel it away from adhesive and proceed to next step.

3. **Non-sterile user:** Apply transducer face directly to adhesive pad.

Sterile user: Grip transducer from below using hand inside cover. Use other hand to deploy cover over transducer and cord. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover. Cover may be repositioned on transducer, as desired.

WARNING


- Ensure no air bubbles are present prior to use. Air left between cover and transducer lens may cause poor image quality.

4. Secure cover with enclosed bands.
5. Inspect cover to ensure there are no holes or tears.

HYDRATING THE COVER

1. **Sterile user:** Use suitable hydration source, such as liquid saline or a saline wipe, to fully saturate patient. To prolong scan time, fully saturate Envision as well.

NOTE: Envision may be hydrated multiple times. Rehydrate as necessary.

2. **Sterile user:** Ensure Envision is fully hydrated prior to use.

STORAGE CONDITIONS

- Avoid storing product in areas of temperature extremes or in direct sunlight.
- Store in a cool, dry place.

DISPOSAL**WARNING**

- *Dispose of single-use components as infectious waste.*

NOTE: For questions or to order additional CIVCO products, please call +1 319-248-6757 or 1-800-445-6741 or visit www.CIVCO.com.

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**Australian
Sponsor Address**

RQSolutions Pty Ltd
44 Del Monte PL.
Copacabana NSW 2251

CIVCO Medical Instruments Co., Inc.
102 First Street South, Kalona, IA 52247-9589 USA
319.248.6757, 800.445.6741
Fax: 319.248.6660, 877.329.2482

COPYRIGHT © 2025 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ENVISION IS A TRADEMARK OF CIVCO MEDICAL SOLUTIONS. 100% GEL-FREE ULTRASOUND PROCEDURES IS A SERVICE MARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.

WWW.CIVCO.COM

2025-01-16 043-859-3J

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MPS GmbH