DISPOSABLE ENDOCAVITY NEEDLE GUIDE

REFERENCE GUIDE

For use with
Canon Medical Systems PVF-651VT, PVM-651VT,
PVN-661VT, PVQ-641V, PVT-661VT;
Siemens EC9-4 (on SONOLINE Antares and S2000)
transducers

REF 810-400

REF 810-400-5

UK
CA
0120

Refer to separate insert for Symbol Glossary.

NEEDLE GUIDE

INTENDED USE

This device provides physicians with a tool for performing needle/instrument guided procedures with the use of the diagnostic ultrasound endocavity transducer.

INDICATIONS FOR USE

- Transvaginal - Diagnostic ultrasound needle/instrument guided procedures such as tissue biopsy, fluid aspiration, catheter placement, treatment, and oocyte retrieval.
- Transrectal - Diagnostic ultrasound needle/instrument guided procedures such as tissue biopsy, fluid aspiration, catheter placement and treatment.

PATIENT POPULATION

- Endocavity Guidance Systems (IVF) - Endocavity Guidance Systems with an IVF indication are for use with transvaginal transducers intended for use on females of all body habitus, including small or slight frame.
- Endocavity Guidance Systems are for use with endocavity (transvaginal or transrectal) transducers intended for use on adults of all body habitus, including small or slight frame.

INTENDED USERS

Endocavity Guidance products should be used by clinicians medically trained in the use and interpretation of ultrasound for interventions. There are no other unique skills or user abilities required for Endocavity Guide use. This may include, but is not limited to: Maternal Fetal physicians, OB/GYN, Reproductive Endocrinologists, Women's Health Nurse Practitioners, Radiologists, Sonographers, Surgeons and Urologists.

PERFORMANCE CHARACTERISTICS

- Endocavity Guidance Systems used during oocyte retrieval must pass MEA and LAL testing to be deemed safe to use in oocyte retrieval procedures.
- Endocavity Guidance Systems are designed to correspond with on-screen system guidelines, directing instruments for accurate placement. Contoured design provides a low-profile, secure fit to the transducer, maximizing patient comfort.

NOTE: For a summary of clinical benefits for this product, visit www.CIVCO.com.

NOTE: One-cell MEA: $\geq 80\%$ embryos developed to expanded blastocyst at 96 hours. USP Endotoxin (LAL) tested and passed with 20 EU or less per device.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

- Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's user guide.
- Use of a needle guide on an endocavity transducer may increase the risk of vaginal or rectal abrasion.
- For illustration purposes only, transducer may be shown without a transducer cover. Always place a cover over transducer to protect patients and user from cross-contamination.
- Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.
- Product is sterilized using ethylene oxide.
- Sterile product is intended for single-use only.
- Do not use if integrity of packaging is violated.
- Do not use if expiration date has passed.
- Do not use if needle guide is damaged or does not fit properly.
- Ensure no air bubbles are present prior to use. Air left between cover and transducer lens may cause poor image quality.
- If the product malfunctions during use or is no longer able to achieve its intended use, stop using the product and call CIVCO.
- Report serious incidents related to the product to CIVCO and the competent authority in your Member State or appropriate regulatory authorities.

NOTE: Product is not made with natural rubber latex.


NEEDLE PATH VERIFICATION

- Before using component for first time, perform needle path verification to verify system and biopsy guide relationships, as described in system user guide.
- See transducer user guide for proper procedure.
- Use care when inserting needle to avoid bending.


NOTE: Use appropriate needle length to reach target area.

USING NEEDLE GUIDANCE


1. Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
2. Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
3. Secure cover with enclosed bands.
4. Inspect cover to ensure there are no holes or tears.


5. Using proper sterile technique, align guide with locating features of transducer.

6. Press down to secure to transducer.

7. If desired, apply secondary cover by following the steps above.

DISPOSAL

WARNING

- *Dispose of single-use components as infectious waste.*
- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*
- *See your system's user guide for reprocessing transducer between use.*

STORAGE CONDITIONS

- Avoid storing product in areas of temperature extremes or in direct sunlight.
- Store in a cool, dry place.

NOTE: For questions or to order additional CIVCO products, please call +1 319-248-6757 or 1-800-445-6741 or visit www.CIVCO.com.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.
Qserve Group UK, Ltd.
RQ Solutions Pty Ltd
CIVCO Medical Instruments Co., Inc.