


TP **PIVOT** PRO™

NEEDLE GUIDE

REFERENCE GUIDE

UK
CA
0120

Refer to separate insert for Symbol Glossary.

NEEDLE GUIDE

INTENDED USE

Device is intended to be used for instrument placement to access anatomical structures under transrectal ultrasound guidance.

COVER

INTENDED USE

Protective cover or sheath placed over an ultrasound transducer / probe / scanhead instruments. The cover is used with the transducer in scanning and needle guided procedures for body surface, endocavity, and intra-operative diagnostic, interventional, or therapeutic ultrasound. The cover also helps to prevent transfer of microorganisms, body fluids, and material to the patient and healthcare worker during reuse of the transducer (both sterile and non-sterile covers). The cover also provides a means for maintenance of a sterile field (sterile covers only).

INDICATIONS FOR USE

- Prostate - Biopsy and minimally invasive puncture.
- Surgical (Prostate)- Biopsy and minimally invasive puncture.

PATIENT POPULATION

TP Pivot Pro is intended for use on adult males of all body habitus, including small or slight frame.

INTENDED USERS

Needle Guides should be used by clinicians medically qualified in the proper use of transperineal ultrasound for guiding needle placement. User groups may include, but are not limited to: Urologist, Urological Surgeon, Urological Nurse, Radiation Oncologist, Physician Assistant and Nurse Practitioner.

PERFORMANCE CHARACTERISTICS

- TP Pivot Pro needle guide provides a secure fit on most transrectal ultrasound transducers with its adjustable locking knob which ensures a stable position allowing users to perform guided, free-handed transperineal prostate biopsies through an initial access needle. TP Pivot Pro offers 5 adjustable heights and varying angles to improve access to different zones throughout the prostate gland. The needle path angle feature decreases the need to angle the probe, thereby reducing probability of image dropout when assessing posterior gland.
- Cover serves as a viral barrier to protect patients, users, and equipment from cross-contamination.

INTENDED CLINICAL BENEFITS

- Freehand transperineal prostate biopsy approach provides detection of clinically significant disease and greater sampling of the anterior prostate.

NOTE: For a summary of clinical claims for this product, visit www.CIVCO.com.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

- Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's user guide.
- Cover and guide are sterilized using ethylene oxide.
- Sterile product is intended for single-use only.
- Do not use if integrity of packaging is violated.
- Do not use if expiration date has passed.
- Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.
- Do not use if needle guide is damaged or does not fit properly.
- Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.
- For illustration purposes only, transducer may be shown without a transducer cover. Always place a cover over transducer to protect patients and user from cross-contamination.
- Ensure no air bubbles are present prior to use. Air left between cover and transducer lens may cause poor image quality.

- Do not overtighten clamping knob. Overtightening may cause transducer damage.
- Ensure introducer needle is visualized when angled, to avoid puncturing the rectal wall.
- Withdraw or remove the introducer needle from the introducer holder to prevent patient injury when moving to a contralateral position.
- If the product malfunctions during use or is no longer able to achieve its intended use, stop using the product and call CIVCO.
- Report serious incidents related to the product to CIVCO and the competent authority in your Member State or appropriate regulatory authorities.
- Use care when handling introducer needle to prevent accidental punctures.

NOTE: Product is not made with natural rubber latex.

NEEDLE PATH VERIFICATION

- Before using device, verify introducer needle is aligned within transducer imaging plane.
- See transducer user guide for proper procedure.
- Use care when inserting needle to avoid bending.

NOTE: Use appropriate needle length to reach target area.


USING NEEDLE GUIDANCE

1. Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
2. Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
3. Secure cover with enclosed bands.
4. Inspect cover to ensure there are no holes or tears.

ATTACHING NEEDLE GUIDE TO TRANSDUCER

NOTE: Guide is compatible with endocavity transducers with a diameter in the range of 16mm - 20mm (0.63 - 0.79 inches).


1. Open transducer clamping knob by turning knob counter-clockwise and pivot knob away from attachment legs. Do not remove knob from guide.
2. Place needle guide over transducer and push down to attach.
3. Align stabilizer window with sagittal image plane transducer lens.
4. Position guide to allow introducer needle to be visible in transducer array once inserted.
5. Pivot transducer clamping knob back into position and turn knob clockwise to securely attach needle guide to transducer. Do not overtighten.

NOTE: If necessary, loosen transducer clamping knob to adjust needle guide position on transducer and ensure needle path is aligned with transducer array. Do not remove knob from guide. Tighten transducer clamping knob to secure needle guide after adjustment.

ADJUSTING NEEDLE GUIDE

1. Determine desired starting needle path height for procedure.
2. Loosen introducer locking knob by turning counter-clockwise approximately one full turn. Do not remove knob from guide tower.
3. Pull introducer holder out of current height position.
4. Move introducer holder up or down to desired needle path height.
5. Push introducer holder into desired needle path height.
6. Tighten introducer locking knob by turning clockwise until introducer holder is securely locked in a horizontal position using alignment detents.


NOTE: Do not remove introducer holder from guide tower.


NOTE: Introducer needle path height and angle adjustments may be made as required during the procedure.

ASSEMBLING INTRODUCER NEEDLE TO NEEDLE GUIDE

1. Using proper sterile technique, remove introducer needle from packaging and remove tip guard from introducer needle.
2. Position needle guide for initial introducer needle placement.
3. Insert introducer needle into introducer holder to desired depth of puncture.

Depth mark indicates approximate insertion distance prior to puncture.


NOTE: Some resistance to introducer needle insertion will be noticed at insertion distance immediately preceding puncture depth. Depth is approximately indicated by depth mark on introducer needle.

NOTE: Introducer needle insertion depth may be adjusted during procedure.

NOTE: Introducer needle accepts 18GA biopsy needles.

ADJUSTING INTRODUCER HOLDER NEEDLE PATH ANGLE

1. Loosen introducer locking knob by turning counter-clockwise approximately one full turn. Do not remove knob from guide tower.
2. Adjust angle of introducer holder to desired angle.
3. If desired, lock angle of needle path by turning introducer locking knob clockwise until secure. Do not overtighten.

**USING NEEDLE GUIDE**

1. Select needle path height and ensure introducer holder is locked into position.
2. Position transducer with needle guide attached against patient's perineal skin for initial prostate biopsy.
3. Insert introducer needle into introducer holder and push needle forward to puncture perineal skin. Position needle to desired depth in subcutaneous tissue.
4. Perform biopsy procedure, adjusting introducer needle path height, angle, and introducer insertion depth as necessary to collect all biopsies required from puncture site.
5. Remove introducer needle from introducer holder.
6. Repeat as required at additional positions for prostate biopsy.

DISASSEMBLING AND REMOVING NEEDLE GUIDE

1. Loosen transducer clamping knob by turning counter-clockwise. Pull guide to remove from transducer.


DISPOSAL**WARNING**

- *Dispose of single-use components as infectious waste.*
- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*
- *See your system's user guide for reprocessing transducer between use.*


STORAGE CONDITIONS

- Avoid storing product in areas of temperature extremes or in direct sunlight.
- Store in a cool, dry place.


NOTE: For questions or to order additional CIVCO products, please call +1 319-248-6757 or 1-800-445-6741 or visit www.CIVCO.com.

- I. Toranj za navođenje
- II. Prozor stabilizatora
- III. Nožice za pričvršćivanje
- IV. Gumb za stezanje sonde
- V. Gumb za zaključavanje uvodnika
- VI. Držač uvodnika

- I. Brava uvodnika sa zaporima za vodoravno poravnanje
- II. Uvodna igla
- III. Platforma za navođenje
- IV. Indikatori visine

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. TP PIVOT PRO IS A TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.
Qserve Group UK Ltd.
RQ Solutions Pty Ltd
CIVCO Medical Instruments Co., Inc.